

HemoStar[®] 
LONG-TERM HEMODIALYSIS CATHETER

HemoStar[®]  **XK**
LONG-TERM HEMODIALYSIS CATHETER

With
Air Guard[®]
VALVED INTRODUCER*



BAIRD
ACCESS SYSTEMS

Advancing the Delivery of Health Care.™

Equipment You Can Count on for...

Exceptional Performance in Long-Term Dialysis and Improved Protection from Blood Loss and Air Embolism



Exceptional Performance

- Polyurethane catheter material – affords strength for longevity and softness for flexibility and patient comfort
- Large lumens and non-restrictive tip design – enable flow rates as high as 500 ml/min⁺⁺

Easy Insertion

- Exceptional kink resistance – simplifies insertion by allowing greater flexibility in tunnel location
- Guidewire channel – on the venous tip facilitates threading for easy over-the-wire insertion
- Malleable tunneler, PTFE introducer and stiff guidewire – facilitate catheter insertion



Optimized Performance for Difficult Patients

HemoStar[®] XK catheter offers the same advantages as the HemoStar[®] catheter with the following features:

- 16 Fr. catheter – allows for flow rates of 500 ml/min⁺⁺ with lower pressures than the 14.5 Fr. HemoStar[®] Catheter
- Increased longevity – larger inner-lumen area allows for flows of 500 ml/min⁺⁺ even when occluded 15%




Improved Protection from Air Embolism and Blood Loss

- 0.15 ml/second air leakage⁺ – most effective valved introducer for preventing air embolism
- Reduction in blood loss – a safer procedure when compared to non-valved introducers

⁺At - 12mm Hg Vacuum Pressure

⁺⁺Blood Pump Setting



HemoStar [®] 14.5 Fr. Catheter Standard Kits				
Product Code	Description	Tip to Cuff	Tip to Bifurcation	
Straight Catheters				
5833150	HemoStar [®] Catheter	15 cm	20 cm	
5833690	HemoStar [®] Catheter	19 cm	24 cm	
5833730	HemoStar [®] Catheter	23 cm	28 cm	
5833270	HemoStar [®] Catheter	27 cm	32 cm	
5833310	HemoStar [®] Catheter	31 cm	36 cm	
5833350	HemoStar [®] Catheter	35 cm	40 cm	
5834420	HemoStar [®] Catheter	42 cm	47 cm	
Alphacurve[®] Catheters				
5835150	HemoStar [®] Catheter	19 cm	25 cm	
5835190	HemoStar [®] Catheter	24 cm	29 cm	
5835230	HemoStar [®] Catheter	28 cm	33 cm	
5835270	HemoStar [®] Catheter	31 cm	37 cm	

HemoStar [®] XK 16 Fr. Catheter Standard Kits				
Product Code	Description	Tip to Cuff	Tip to Bifurcation	
Straight Catheters				
5883690	HemoStar [®] XK Catheter	19 cm	24 cm	
5883730	HemoStar [®] XK Catheter	23 cm	28 cm	
5883270	HemoStar [®] XK Catheter	27 cm	32 cm	
5883310	HemoStar [®] XK Catheter	31 cm	36 cm	
5883350	HemoStar [®] XK Catheter	35 cm	40 cm	
5884420	HemoStar [®] XK Catheter	42 cm	47 cm	
Alphacurve[®] Catheters				
5885150	HemoStar [®] XK Catheter	19 cm	25 cm	
5885190	HemoStar [®] XK Catheter	24 cm	29 cm	
5885230	HemoStar [®] XK Catheter	28 cm	33 cm	
5885270	HemoStar [®] XK Catheter	31 cm	37 cm	

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

x _____

Physician's Signature



ACCESS SYSTEMS

Bard Access Systems, Inc.

605 North 5600 West
Salt Lake City, Utah 84116 USA
801-522-5000
www.bardaccess.com

Customer Service: 800-545-0890
Clinical Information: 800-443-3385